

# Tidy Tubing™

Retractable Oxygen Tubing for  
Fall Prevention & Convenience

Model 1834



[CLICK TO LEARN MORE](#)



## Extends to 10 Feet

### RECOMMENDED USE:

- Portable oxygen concentrator(POC)
- Pull Cart
- Oxygen Tank Bag
- Wheelchair



### WHY TIDY TUBING?

- Eliminates the need for long O2 tubing that can get hung up, tangled, or kinked.
- Fantastic for fall prevention

### INCLUDED:

- Clothing clip to prevent nose tug
- 2 foot Adult Vanish Cannula
  - soft&pliable, reduces irritation
- Swivel connectors on each end to prevent tangling



### ADDITIONAL MODELS:

- 1833: extends to 5 feet
- 1835: extends to 15 feet

